1 Intended use

Intended purpose

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Patient target groups

- Patients with permanent teeth
- Patients with deciduous teeth

Intended users / Special training

- Dentists
- No special training required.

Use

For dental use only.

Description

Cention® N is a self-curing, radiopaque ($\geq 200\%$ Al) filling material with light-curing option, which releases fluoride, calcium and hydroxide ions (according to ISO 4049:2019 Type 1, Class 3). It is suitable for use as a full volume replacement material in the restoration of Class I, II and V cavities. Cention N is also suitable for restoring occlusal surfaces.

100% aluminum has a radiopacity equivalent to that of dentin and 200% aluminum is equivalent to enamel.

Cention N can be optionally cured with light in the wavelength range of 400 – 500 nm. Cention N is intended to be manually mixed.

Indications

Missing tooth structure in anterior and posterior teeth

Types of restorations:

- Suitable for Class I (occlusal) and II (occlusal and proximal) restorations in conjunction with an adhesive.
- Suitable for use without an adhesive in Class I and Class II cavities with retentive preparation, provided that the distance between the cavity margins and the cusp tip is ≥ 1 mm.
- Class V restorations (cervical) in permanent teeth in conjunction with an adhesive.

Contraindications

The use of the product is contraindicated if the patient is known to be allergic to any of its ingredients.

Limitations of use

- Due to esthetic reasons, Cention N is not suitable for Class III and IV restorations.
- If a dry working field cannot be established.
- If the stipulated working technique cannot be applied.
- When used without an adhesive to merely fill the proximal box (slot cavity).
- When used without an adhesive to replace a missing cusp.
- When used as a luting composite.

Side effects

In rare cases, ingredients of Cention N may cause a sensitizing reaction. The product must not be used in such cases. In order to avoid irritation of the pulp, provide areas close to the pulp with suitable pulp/dentin protection. Selectively apply a calcium hydroxide-based preparation in areas close to the pulp and cover with a suitable cavity liner, if necessary.

Interactions

Substances containing eugenol / clove oil may inhibit the polymerization of methacrylate-based materials. Consequently, the application of such materials together with Cention N must be avoided. Discolouration may occur in combination with cationic mouthwashes, plaque disclosing agents and chlorhexidine. Disinfectants with an oxidative effect (e.g. hydrogen peroxide) may interact with the initiator system, which in turn may impair the curing process.

Clinical benefit

Reconstruction of chewing function

Composition

Cention N Liquid contains UDMA, aromatic-aliphatic UDMA, DCP, PEG-DMA, initiators and mint flavour.

Cention N Powder contains calcium fluoro-silicate glass, copolymer, barium-aluminium silicate glass, calcium-barium-aluminium fluoro-silicate glass, barium glass, and ytterbium trifluoride.

Once mixed (Powder:Liquid = 4.6:1 parts by weight, which corresponds to 1 measuring scoop of Powder + 1 drop of Liquid), Cention N contains 58-59 vol% of inorganic fillers.

Particle size of the inorganic fillers: between 0.15 and 1.7 μm .

2 Application

I. Isolation

Adequate relative or absolute isolation is required.

II. Cavity preparation

Cavity preparation when using Cention N without an adhesive:

- When using Cention N without an adhesive, the preparation guidelines for the placement of amalgam fillings apply (retentive preparation design, resistance form). Retentive undercuts need to be prepared in all areas of the cavity. Enamel portions not supported by dentin must be avoided. In order to prevent fractures of class II fillings, the area of the isthmus should be sufficiently deep and wide: on the occlusal/proximal aspect it should be at least 1.5 mm in depth, and the isthmus width should be at least 1.5 mm. The cavity floor must be rounded from the occlusal to the proximal aspect in order to reduce tension peaks. A dental adhesive should always be used when large defects are restored. Preparations with severely undermined areas that lead to a weakening of the cusps should be avoided. The distance between the cavity margins and the cusp tips should be ≥ 1 mm. Enamel margins should not be bevelled.

Cention® N

[en]	Instructions for Use Polymer-based dental restorative material (intraoral dual-curing)	[pt]	Instruções de Uso Material restaurador odontológico baseado em polímero (polimerização dual intraoral)
[fr]	Mode d'emploi Matériau de restauration dentaire à base de polymères (polymérisation duale intraorale)	[tr]	Kullanma Talimatı Polimer bazlı dental restoratif malzeme (ağız içinde çift sertleşen)
[es]	Instrucciones de uso Material dental restaurador basado en polímeros (polimerización dual intraoral)	[ru]	Инструкция по применению Полимерный стоматологический пломбирочный материал (двойное отверждение в полости рта)

Subsequently, rinse the cavity with water spray to remove all the residue and dry with water- and oil-free air.

Cavity preparation when optionally using Cention N with an adhesive:

- The cavity is prepared according to the principles of the adhesive technique, i.e. by preserving as much of the tooth structure as possible. Do not prepare sharp internal angles or additional undercuts in caries-free areas. The geometry of the cavity is predetermined by the dimensions of the caries lesion or the old filling. Slightly bevel or round out the enamel margins using finishing diamonds (grit size 25 – 40 µm). Subsequently, rinse the cavity with water to remove all residue and dry with water- and oil-free air. Note: Do not overdry the dentin.

III. Pulp protection / Base

- If no adhesive is used, areas close to the pulp should be selectively coated with a calcium hydroxide preparation.
- If a dental adhesive is used, a base material does not have to be applied. Only in very deep cavities, selectively coat the area close to the pulp with a calcium hydroxide-based preparation.

IV. Placement of matrix / interdental wedge

Use a wrap-around matrix for cavities affecting the proximal area or a sectional matrix band and secure it with wedges.

V. Conditioning / Applying the adhesive

When using Cention N without an adhesive, do not condition the tooth structure with phosphoric acid. Condition and apply the adhesive according to the Instructions for Use of the product in use.

VI. Dispensing and mixing

- The mixing ratio for a small cavity is:
1 measuring scoop of Powder and 1 drop of Liquid (this corresponds to a weight ratio of 4.6:1).
Note: use the supplied measuring spoon.
- For larger cavities, a mixing ratio of either 2 measuring scoops of Powder and 2 drops of Liquid or 3 measuring scoops of Powder and 3 drops of Liquid should be used.
- Shake the Powder bottle well before use.
- Dispense Powder and Liquid next to each other on a mixing pad.
- Hold the bottle containing the liquid perpendicular when dispensing and squeeze to extrude one drop. It is important to make sure that the drop detaches before it makes contact with the mixing pad and to avoid the inclusion of air bubbles.
- Overfill the measuring scoop with Powder and scrape off the excess material on the plastic insert in the bottle neck. Close the Powder bottle and the Liquid bottle carefully after use.
- Separate the powder into two equally large portions using a plastic spatula. Spread the Liquid to enlarge the surface.
- Mix the first portion of Powder with the entire Liquid dispensed on the mixing pad. Once the components have been thoroughly mixed, add the remaining Powder and mix again until a homogeneous consistency is achieved (45 – 60 s). The working time is 3 min from the start of mixing.
- Apply the material to the cavity using a suitable instrument. Carefully adapt and condense it and remove any occlusal excess.

Processing and setting times

Mixing time (at 23 °C)	– 45 – 60 seconds
Working time	--- 3 minutes
Setting time (self-curing)	----- 5 minutes (from the start of mixing)

The indicated working and setting times were determined under clinical conditions. High temperatures during mixing decrease the working and setting time; low temperatures increase it.

- When using the self-curing mode, the matrix band can be removed after the setting time has elapsed and finishing can be started.
- In order to remove the matrix band and start finishing of the restoration in a timely fashion, the restoration can be optionally light-cured after placement. For the recommendations regarding exposure time (Exposure time) per increment and light intensity (Light intensity) see table 1 (Table 1).

When filling cavities measuring more than 4 mm in depth, wait until the entire setting time has elapsed.

VII. Finishing / Checking the occlusion

Remove any excess material with tungsten carbide or diamond finishers after polymerization. Check the occlusion and articulation and apply appropriate grinding corrections to prevent premature contacts or undesired articulation paths on the surface of the restoration. If needed, use polishers as well as polishing discs and polishing strips to polish the restoration to a high gloss.

Additional notes

- In the case of adjustments, additional Cention N can be directly applied to the polymerized material. If the inhibition layer has already been removed, the restoration surface must be roughened and coated with adhesive before fresh Cention N can be applied. The adhesive must be applied according to the manufacturer's instructions.
- Use Cention N at room temperature. Cold components may be difficult to mix.
- Do not use oxidizing disinfectants to disinfect the bottles.

3 Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclar.com, and your responsible competent authority.
- The current Instructions for Use and the explanation of symbols are available on the website: www.ivoclar.com/eIFU
- The Summary of Safety and Clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>.
- Basic UDI-DI: 76152082AFILL016JP

Warnings

- Observe the Safety Data Sheet (SDS) (available at Ivoclar Vivadent AG www.ivoclar.com).
- Unpolymerized Cention N should not come in contact with skin / mucous membrane or eyes. Unpolymerized Cention N may have a slight irritating effect and may lead to a sensitization against methacrylates. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Disposal information

Remaining stocks or removed restorations must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks. The following known residual clinical risks exist:

- Failure of the adhesive bond (loss of the filling)
- Inclusion of air bubbles during the placement of the filling
- Chipping, fractures
- Ingestion of material
- Postoperative sensitivity

4 Shelf life and storage

- Storage temperature 2–28 °C (36–82 °F)
If the storage temperature is exceeded, the proper function of the product can be impaired.
- Do not use the product after the indicated date of expiration.
- Expiry date: see information on the Liquid and Powder bottles and packages.
- Once the Cention N Powder bottle has been opened, the product should be used up within 8 months.
- Close bottles immediately after use; exposure to light causes premature polymerization.
- Before use, visually inspect the packaging and the product for damage. In case of any doubts, please contact Ivoclar Vivadent AG or your local dealer.

5 Additional information

Keep material out of the reach of children!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

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