Monobond® N

[en] Instructions for Use

Universal primer

Rx ONLY

CE 0123

Date information prepared:
2023-11-23 / Rev. 0

Manufacturer:
Ivoclar Vivadent AG
Bendererstrasse 2
9494 Schaan/Liechtenstein
www.ivoclar.com

ivoclar

English

1 Intended use

Intended purpose

Preparing indirect dental restorations for adhesive cementation

Patient target group

- Patients with permanent teeth
- Adult patients with dental implants

Intended users

- Dentists (clinical procedure)
- Dental assistants (preparation and extraoral application)

Special training

No special training required.

Use

For dental use only.

Description

Monobond® N is a universal primer promoting an adhesive bond between luting composites (particularly the Variolink and Multilink product lines) and indirect restorative materials.

Indications

- Missing tooth structure in anterior and posterior teeth
- Partial edentulism in the anterior and posterior region

Areas of application/types of restorations

The product serves as bonding agent and is used to create a durable chemical bond between luting composites and glass/oxide ceramic, metal and fibre-reinforced composite restorations.

Contraindications

The use of the product is contraindicated if the patient is known to be allergic to any of its ingredients.

Limitations of use

- Do not use if a dry working field cannot be established or the prescribed application technique cannot be applied.
- Contamination with saliva, blood or water must absolutely be avoided during the entire procedure. In case of contamination, the entire treatment of the restoration must be repeated.
- In order to achieve an optimum bond, do not clean zirconium oxide and metal surfaces with phosphoric acid.
- Do not insert the brush into the Monobond N bottle. This may lead to contamination.
- Carefully close the Monobond N bottle immediately after use to prevent evaporation of the solvent.

Side effects

In rare cases, components of Monobond N may lead to sensitization. In these cases, Monobond N should no longer be used.

Interactions

There are no known interactions to date.

Clinical benefit

- Reconstruction of chewing function
- Restoration of esthetics

Composition

Alcohol solution of silane methacrylate, phosphoric acid methacrylate and sulphide methacrylate.

2 Application

For more detailed information, also refer to the Instructions for Use of the products used in conjunction with Monobond N.

Preparation of the restoration

a) Glass-ceramic restorations (e.g. IPS Empress®)

- Etch the restoration with 5% hydrofluoric acid (e.g. IPS® Ceramic etching gel) for 60 s or according to the instructions of the manufacturer of the restorative material.
- Thoroughly rinse the restoration with water spray and dry with oil-free air.

b) Lithium disilicate glass-ceramic restorations (e.g. IPS e.max® Press, IPS e.max CAD)

- Etch the restoration with 5% hydrofluoric acid (e.g. IPS Ceramic etching gel) for 20 s or according to the instructions of the manufacturer of the restorative material.
- Thoroughly rinse the restoration with water spray and dry with oil-free air.

c) Zirconium oxide restorations (e.g. IPS e.max ZirCAD) or aluminium oxide ceramic restorations

- Sandblast the bonding surface (sandblasting parameters according to the instructions of the manufacturer of the restorative material).
 - If necessary, clean the restoration in an ultrasonic unit for about one minute.
 - Thoroughly rinse the restoration with water spray and dry with oil-free air.
- IMPORTANT!** In order to create a strong bond, do not clean the zirconium oxide surfaces with phosphoric acid.

d) Metal or metal-supported restorations

- Sandblast the restoration surface (sandblasting parameters according to the instructions of the manufacturer of the restorative material) until an even, matt or non-glossy surface has been achieved.
 - If necessary, clean the restoration in an ultrasonic unit for about one minute.
 - Thoroughly rinse the restoration with water spray and dry with oil-free air.
- IMPORTANT!** In order to create a strong bond, do not clean the metal surfaces with phosphoric acid.

e) Fibre-reinforced composite restorations

- Condition the restoration according to the instructions of the manufacturer of the restorative material.

After having prepared the restoration, proceed as follows:

- Apply a thin coat of Monobond N with a brush or a microbrush to the pre-treated surfaces (avoid pooling when treating crowns). Allow the material to react for 60 seconds. Subsequently, disperse any remaining excess with a strong stream of air.
- Apply the luting composite according to the respective Instructions for Use and insert the restoration.

3 Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclar.com and your responsible competent authority.
- The current Instructions for Use and the explanation of symbols are available on the website: www.ivoclar.com/eIFU.
- The Summary of Safety and Clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>.
- Basic UDI-DI: 76152082ABOND009J2

Warnings

- Observe the Safety Data Sheet (SDS) (available at www.ivoclar.com).
- The product contains ethanol. Keep material away from sources of ignition.
- Prevent any contact of uncured material with the skin/mucous membrane and eyes.
- Uncured material may have a slight irritating effect and may lead to a sensitization against methacrylates.
- Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Disposal information

Remaining stocks or removed restorations must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks.

The following clinical residual risks are known:

- Failure of the adhesive bond

4 Shelf life and storage

- Storage temperature: 2-28 °C (36-82 °F)
- Do not use the product after the indicated date of expiration.
- Expiry date: See note on bottle or packaging.
- Before use, visually inspect the packaging and the product for damage. In case of any doubts, please contact Ivoclar Vivadent AG or your local dealer.

5 Additional information

Keep material out of the reach of children!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the product for its suitability and use for any purpose not explicitly stated in the Instructions.

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