Variolink® N

[en] Instructions for Use

Polymer-based dental luting material (intraoral dual-curing / intraoral light-curing)

C E 0123

Manufacturer:
Ivoclar Vivadent AG
Bendererstrasse 2
9494 Schaan/Liechtenstein
www.ivoclar.com

Date information prepared:
2024-02-06/Rev. 2
web

ivoclar

English

1 Intended Use

Intended purpose

Adhesive cementation of indirect anterior and posterior restorations

Patient target groups

Patients with permanent teeth

Intended users / Special training

- Dentists
- No special training required

Use

For dental use only.

Description

Variolink N is a radiopaque, dual-/light-curing luting composite system for the adhesive luting of glass-ceramic, lithium disilicate glass-ceramic and composite restorations (according to ISO 4049:2019 Type 2, Class 3). The special filler composition gives Variolink N a very high radiopacity ($\geq 300\%$ Al). 100% aluminum has a radiopacity equivalent to that of dentin and 200% aluminum is equivalent to enamel. Variolink N can be intraorally cured with light in the wavelength range of 400–500 nm.

Shades

Variolink N is available in 4 base shades and 2 catalyst shades as well as 2 levels of viscosity:

Base	Shade A–D
Bleach XL	BL 1
Transparent	T
White	A1
Yellow (universal)	A3
Catalyst high/low viscosity	Shade A–D
Transparent	T
Yellow (universal)	A3

We recommend the base with the low-viscosity catalyst paste for the placement of crowns and partial crowns, the base in combination with the high-viscosity catalyst paste for inlays/onlays and veneers. The dentist, however, is free to choose the proper viscosity for each indication.

In order to evaluate the overall effect of the restoration in conjunction with the various Variolink N shades prior to permanent cementation, we recommend using the Variolink N Try-In Pastes.

Working time

The working and curing times depend on the ambient temperature. The following times apply as soon as Variolink N has been started to be hand-mixed:

When used purely self-curing	Room temperature 23°C ± 1°C/ 73°F ± 2°F	Intraoral 37°C ± 1°C/ 99°F ± 2°F
Working Time	5 min.	3.5 min.
Curing time (incl. working time)	12 min.	7 min.

Mixing ratio

Mix base paste and catalyst in a 1:1 ratio.

Indications

- Missing tooth structure in anterior and posterior teeth
- Partial edentulism in the anterior and posterior region

Areas of application

Variolink N is a dual-curing luting composite system for the adhesive luting of glass-ceramic, lithium disilicate glass-ceramic and composite restorations (inlays, onlays, partial crowns, crowns, bridges, veneers and endodontic posts).

Contraindications

The use of this product is contraindicated if the patient is known to be allergic to any of its ingredients.

Limitations of use

Do not use if a dry working field cannot be established or the prescribed application technique cannot be applied.

Side effects

Systemic side effects are not known to date.

Interactions

- Phenolic substances (e.g. eugenol, wintergreen oil) inhibit polymerization. Consequently, the application of materials containing these substances, e.g. mouth rinses and temporary cements, is to be avoided.
- Disinfectants with an oxidative effect (e.g. hydrogen peroxide) may interact with the initiator system, which in turn may impair the curing process. Therefore, do not disinfect the preparation using oxidative agents.
- Alkaline jet media applied on dentin can compromise the effect of self-etching adhesives.
- Hemostatic substances may inhibit polymerization and/or lead to discoloration. Therefore, the directions for use of these substances must be observed.

Clinical benefit

- Reconstruction of chewing function
- Restoration of esthetics

Composition

Variolink N is composed of barium glass, ytterbiumtrifluoride, Bis-GMA, UDMA, TEGDMA, Si-Zr mixed oxide, barium-aluminium-fluorosilicate glass.

Filler particle sizes range from 0.15 to 15.5 µm.

Total content of inorganic filler (percentage by volume):

Base	Filler %vol.
	46.7
Catalyst high/low viscosity	
low viscosity	43.6
high viscosity	52.0

2 Application

For more detailed information also refer to the instructions for use of the products used in conjunction with Variolink N.

2.1 Removal of the temporary restoration and cleaning of the cavity

Remove any possible residue of the temporary luting composite from the cavity or the prepared tooth with a polishing brush and an oil- and fluoride-free cleaning paste. Rinse with water spray. Subsequently, dry with oil- and moisture-free air. Avoid overdrying.

Note: Cleaning with alcohol can lead to dehydration of dentin.

2.2 Try-in of the restoration and isolation

Seat the restoration using the desired Try-In Paste (e.g. Variolink N Try-In Paste) and check the shade, fit and occlusion of the restoration. Care should be taken when checking the occlusion of fragile and brittle ceramic objects before they are permanently cemented, as there is a risk of fracture. If necessary, make adjustments with fine diamonds at medium speed and with light pressure and adequate water cooling. Polish ground surfaces. Adequate relative or absolute isolation using suitable auxiliaries, such as OptraGate® or OptraDam Plus, is required. If the dental hard tissue is contaminated with blood or saliva during the try-in procedure, it needs to be cleaned again as described in section 2.1.

2.3 Pre-treatment of the restoration surfaces

- 2.3.1 If a restoration which has been pre-treated in the dental laboratory is tried in, the bonding surface of the restoration needs to be prepared as follows after the try-in:
- Following try-in, thoroughly rinse the restoration with water spray and dry with oil-free air.
 - Shake the cleaning agent (e.g. Ivoclean) before use and completely coat the bonding surface of the restoration with the product using a microbrush or brush.
 - Leave the cleaning agent (e.g. Ivoclean) to react for 20 s, then thoroughly rinse with water spray and dry with oil-free air.
 - Next, prime the bonding surface of the restoration with a suitable bonding agent (e.g. Monobond® N). Make sure to observe the Instructions for Use of the bonding agent used.
- 2.3.2 Restorations that have not been pre-treated by the dental laboratory should be rinsed with water spray and dried after the try-in procedure.
- 2.3.2.1 Subsequently, condition the bonding surface of the restoration as follows:
- a) Glass-ceramic restorations (e.g. IPS Empress®)
 - Etch the restoration with 5% hydrofluoric acid (e.g. IPS® Ceramic etching gel) for 60 s or according to the instructions of the manufacturer of the restorative material.
 - Thoroughly rinse the restoration with water spray and dry with oil-free air.
 - b) Lithium disilicate glass-ceramic restorations (e.g. IPS e.max® Press, IPS e.max CAD)
 - Etch the restoration with 5% hydrofluoric acid (e.g. IPS Ceramic etching gel) for 20 s or according to the instructions of the manufacturer of the restorative material.
 - Rinse the restoration thoroughly with water spray and dry with oil-free air.
 - c) Composite or fibre-reinforced composite restorations:
 - Blast the restoration surfaces (according to the instructions of the manufacturer of the restorative material).
 - If necessary, clean the restoration in an ultrasonic unit for about 1 minute.
 - Thoroughly rinse with water spray and dry with oil-free air.
- 2.3.2.2 Apply primer (e.g. Monobond N) to the pre-treated surfaces using a brush or microbrush, let it react for 60 s and then disperse with a strong stream of air.

Important! Contamination with saliva, blood or water must absolutely be avoided during the entire procedure. In case of contamination, the surfaces of glass-ceramic, lithium disilicate glass-ceramic and composite restorations can be etched with phosphoric acid. Subsequently, the primer needs to be applied again.

2.4 Pre-treatment of the preparation and application of the adhesive

Please observe the instructions for use of the adhesive used.

- 2.4.1 Conditioning with phosphoric acid gel (optional): The bond to enamel can be improved by selectively etching the enamel or by applying the etch & rinse technique. Unprepared enamel surfaces must be conditioned with phosphoric acid gel (e.g. N-Etch). Please observe the instructions for use of the phosphoric acid gel.
- a) Selective enamel etching
 - Apply the phosphoric acid gel onto the enamel and allow it to react for 15–30 s. Then rinse thoroughly with a vigorous stream of water for at least 5 s and dry with compressed air until the etched enamel surfaces appear chalky white.
 - b) Etch & rinse technique
 - Apply the phosphoric acid gel onto the prepared enamel first, and then to the dentin. The etchant should be left to react on the enamel for 15–30 s and on the dentin for 10–15 s. Then rinse thoroughly with a vigorous stream of water for at least 5 s and dry with compressed air until the etched enamel surfaces appear chalky white.
- 2.4.2 Application of the adhesive (e.g. Tetric N-Bond Universal)
- Starting with the enamel, completely coat the tooth surfaces to be treated with Tetric N-Bond Universal.
 - The adhesive must be scrubbed into the tooth surface for at least 20 s. This time must not be shortened. Applying the adhesive on the tooth surface without scrubbing is inadequate.
 - Disperse the adhesive with oil- and moisture-free compressed air until a glossy, immobile film layer results.
- Important information:**
- Avoid pooling, since this can compromise the fitting accuracy of the permanent restoration.
 - Exposure to intensive light should be avoided during application.
 - Prevent any contamination of the applicator/cannula with blood, saliva, sulcus fluid or water during application. In case of contamination, the application must be repeated with a new applicator/new cannula.
- 2.4.3 Light-curing the adhesive
- Tetric N-Bond Universal can be light-cured for 10 s at a light intensity of 500 to 1,400 mW/cm².
- 2.5 **Mixing of Variolink N Base and Catalyst**
- Mix Variolink N Base and Catalyst in a 1:1 ratio on a mixing pad for 10 s immediately before application.
 - Too intensive mixing can cause air to become trapped in the paste. This may delay polymerization and compromise mechanical properties.
 - Variolink N is a light-/dual-curing material and thus sensitive to blue light (operatory light, ambient light). Avoid intensive light during application.
 - The working time of the mixed Variolink N is about 3.5 min at a temperature of 37 °C/99 °F.
- 2.6 **Seating of the restoration**
- 2.6.1 Inlays, onlays, crowns and bridges
- Apply the mixed Variolink N to the area of the restoration which is to be bonded.
 - First, place restoration in situ with slight pressure and remove excess with a brush or another suitable instrument.
 - When working on large restorations, make sure to remove excess material in time, especially in areas that are difficult to reach (proximal or gingival margins).
 - Increase pressure and maintain for several seconds. Remove the Variolink N excess thus forced out with a brush.
- 2.6.2 Veneers
- Apply Variolink N Base with the application tip or a brush/spatula directly to the internal surface of the restoration.
- Note:** Direct application of Variolink N Base to a preparation that has been pre-treated with adhesive may – depending on the adhesive used – lead to an acceleration of the curing process and affect the accuracy of fit of the restoration.
- Seat the restoration exerting light pressure and remove gross excess with a suitable instrument (e.g. spatula, brush).
 - Make sure to remove excess material within the processing time, especially in areas that are difficult to reach (proximal areas and gingival margins).
 - Increase pressure and hold for a few seconds.

- Maintain the pressure and light-cure the central area of the veneer for 3–4 s at a light intensity of approx. 650 mW/cm² (e.g. Bluephase C8, LOW mode) to tack the restoration into place.
Important: Do not light-cure any proximal or marginal regions directly!
- Remove further excess with a suitable instrument. When doing this, take care not to wipe any material out of the marginal areas.

2.7 Polymerization

- As with all composite systems, Variolink N is subject to oxygen inhibition. In other words, the surface layer (approx. 50 µm) does not polymerize during the curing process, as it comes in contact with atmospheric oxygen. To prevent this, cover the restoration margins with glycerine gel/airblock (e.g. Liquid Strip) immediately after excess removal.
- Polymerize Variolink N segment by segment. Begin with the proximal margins.

Light intensity	Curing time per mm ceramic and segment
500–800 mW/cm ²	20 s
900–1300 mW/cm ²	10 s

- Rinse off glycerine gel/airblock (e.g. Liquid Strip)

2.8 Finishing of the restoration

- Check occlusion and functional movements and make adjustments if necessary.
- Finish the cement lines with finishing diamonds if necessary.
- Smooth out the cement lines using finishing and polishing strips and polish them with suitable polishing instruments (e.g. OptraGloss®).
- If necessary, finish the restoration margins with suitable polishers (e.g. OptraGloss).

3 Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclar.com, and your responsible competent authority.
- The current Instructions for Use are available in the download section of the Ivoclar Vivadent AG website (www.ivoclar.com).
- Explanation of symbols: www.ivoclar.com/eIFU
- The Summary of Safety and Clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>.
- Basic UDI-DI: 76152082ACEME001DV

Warnings

- Observe the Safety Data Sheet (SDS) (available in the download section of the Ivoclar Vivadent AG website www.ivoclar.com).
- Prevent any contact of uncured Variolink N with the skin/mucous membrane and eyes.
- Uncured Variolink N may cause slight irritation and may lead to a sensitization against methacrylates.
- Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Disposal information

Remaining stocks must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks. The following known residual clinical risk exists:

- Failure of the adhesive bond

4 Shelf life and storage

- Storage temperature
 - Variolink N Base: 2–28 °C/36–82 °F.
 - Variolink N Catalyst: 2–8 °C/36–46 °F (cool storage)
- Close Variolink N syringes immediately after use. Exposure to light causes premature polymerization.
- Do not use the product after the indicated date of expiration.
- Date of expiration: see information on syringes and packages
- Before use, visually inspect the packaging and the product for damage. In case of any doubt, consult Ivoclar Vivadent AG or your local trade partner.

5 Additional information

Keep material out of the reach of children!

Not all products are available in all countries.

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

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