Multilink® Hybrid Abutment

[en] Instructions for Use

Polymer-based dental luting material (extraoral self-curing)

Rx ONLY

C € 0123

Date information prepared:
2023-07-19 / Rev. 0


Manufacturer:
Ivoclar Vivadent AG
Bendererstrasse 2
9494 Schaan/Liechtenstein
www.ivoclar.com


ivoclar

English

Hybrid abutment

Hybrid abutment crown

1 Intended use

Intended purpose

Permanent extraoral bonding of structures made from lithium disilicate glass-ceramic, zirconium oxide or PMMA on bonding bases made from titanium/titanium alloy for the fabrication of implant-supported hybrid restorations for the replacement of single teeth.

Patient target group

Adult patients with dental implants

Intended users / Special training

- Dentists (chairside fabrication of restorations)
- Dental laboratory technicians (fabrication of restorations in the dental laboratory)
- No special training required.

Use

For dental use only.

Description

Multilink Hybrid Abutment is a polymer-based dental luting composite (extraoral self-curing) designed for the permanent cementation of structures made from lithium disilicate glass-ceramic, zirconium oxide or PMMA on bonding bases made from titanium/titanium alloy for the fabrication of hybrid abutments or hybrid abutment crowns.

Technical specifications

Classification acc. to ISO 4049

- Type 2 / Class 1
- Polymer-based dental luting material (extraoral self-curing)

Shades

HO 0 (High Opacity)

Working time

The working and curing times depend on the ambient temperature. Once Multilink Hybrid Abutment has been dispensed from the automix syringe, the following times apply:

	At room temperature 23 °C ± 1 °C (73 °F ± 2 °F)
Working time	approx. 2 min
Curing time (including working time)	approx. 7 min

The working and curing times become shorter at higher temperatures, while they become longer at lower temperatures.

Mixing ratio

Multilink Hybrid Abutment is always dispensed from the automix syringe in the optimum ratio.

Indications

Partial edentulism in the anterior and posterior region

Areas of application:

Permanent extraoral bonding of structures made from lithium disilicate glass-ceramic, zirconium oxide or PMMA to bonding bases made from titanium/titanium alloy for the fabrication of implant supported hybrid restorations for the replacement of single teeth.

Contraindications

The product is contraindicated

- for intraoral luting;
- if the stipulated working technique cannot be applied;
- if the patient is known to be allergic to components of Multilink Hybrid Abutment.

Limitations of use

Multilink Hybrid Abutment should be processed at room temperature. Low(er) temperatures may cause difficulty in dispensing and mixing the material and may prolong the working and curing time.

Side effects

Systemic side effects are not known. In individual cases, allergic reactions to individual components have been reported.

Interactions

- Phenolic substances (e.g. eugenol, wintergreen oil) inhibit polymerization. Consequently, the application of products containing these components, e.g. mouth rinses and temporary cements, must be avoided.
- Disinfectants with an oxidative effect (e.g. hydrogen peroxide) may interact with the initiator system, which in turn may impair the curing process. Therefore, do not disinfect the preparation using oxidative agents.

Clinical benefit

- Reconstruction of chewing function
- Restoration of esthetics

Composition

Barium glass, Bis-EMA, titanium oxide, ytterbium trifluoride, HEMA, Bis-GMA, UDMA, Si-Zr mixed oxide, barium aluminium fluorosilicate glass, silicon dioxide

The total content of inorganic filler is approx. 36 vol%.
Particle size of inorganic fillers: between 0.15 µm and 7.6 µm

2 Usage

For more detailed information also refer to the Instructions for Use of the products used in conjunction with Multilink Hybrid Abutment.

2.1 Preparation of the bonding base (e.g. Viteo® Base Ti)

The bonding base should be prepared according to the instructions of the manufacturer.

- The bonding base is cleaned in an ultrasonic bath or with a steam cleaner and then dried with blown air.
- The bonding base is screwed to the model analog.
- The structure is placed on the bonding base and the relative position is marked with a waterproof pen. This makes it easier to attain the correct relative end position when the parts are assembled at a later stage.
- The emergence profile of the bonding base must not be blasted or modified in any way!
- **If the manufacturer recommends that the contact surfaces be sandblasted, observe the following procedure:**
 - Apply hard modelling wax to protect the emergence profile, since this type of material is easy to remove later on.
 - Also seal the screw channel with wax.
 - Carefully blast the bonding area surface according to the instructions of the manufacturer.
 - Use an instrument and a steam jet for cleaning. Make sure to remove any wax residue meticulously.
- Clean the base in an ultrasonic bath or with a steam cleaner and then dry it with blown air.
- After the bonding surface has been cleaned, it must not be contaminated under any circumstances as this would impair the bond.
- Apply a suitable bonding agent for titanium (e.g. Monobond Plus) on the cleaned bonding surfaces and allow to react for 60 s. After the reaction time, dry the remaining residue with water- and oil-free air.
- Seal the screw channel with a foam pellet or wax.
- The bonding surface must not be contaminated in the process.
- If using a Viteo Base Ti, insert the Viteo Screw Channel Pin into the screw channel of the titanium base. The Channel Pin can be shortened with a scalpel.

2.2 Preparation of the structures

2.2.1 made from lithium disilicate glass-ceramics (e.g. IPS e.max® CAD/Press)

- The ceramic structure must not be blasted when it is prepared for cementation.
- Clean the ceramic structure in an ultrasonic bath or with a steam cleaner and then dry it with blown air.
- After the bonding surface has been cleaned, it must not be contaminated under any circumstances as this would impair the bond.
- Wax can be applied to protect the outer surfaces and the glazed areas.
- There are two options to prepare the bonding surface:
 - a) Etch the bonding surface with 5% hydrofluoric acid (IPS Ceramic etching gel) for 20 s. Subsequently, thoroughly rinse the bonding surface under running water and then dry it with oil-free air. Apply a silane-containing bonding agent (e.g. Monobond Plus) to the cleaned bonding surface and allow to react for 60 s. After the reaction time, dry the remaining residue with water- and oil-free air.
 - b) Apply Monobond Etch & Prime® on the bonding surface using a microbrush and scrub it into the surface for 20 s. Allow to react for another 40 s. Then thoroughly rinse off Monobond Etch & Prime with water and dry the restoration with a strong stream of water- and oil free air for approximately 10 seconds.

2.2.2 made from zirconium oxide (ZrO₂)

- In order to facilitate its handling, attach the ceramic structure to a holder (e.g. microbrush) before blasting it.
- Protect the outer surfaces of the ceramic structure with a coating of wax.
- In order to better control the blasting procedure, colour the bonding surface with a permanent marker.
- Carefully blast the bonding surface of the ceramic structure according to the instructions of the manufacturer.
- Clean the ceramic structure with a steam cleaner or in an ultrasonic bath.
- After the bonding surface has been cleaned, it must not be contaminated under any circumstances as this would impair the bond.
- Apply a suitable bonding agent for zirconium oxide (e.g. Monobond® Plus) on the cleaned bonding surfaces and allow to react for 60 s. After the reaction time, dry the remaining residue with water- and oil-free air.

2.2.3 made from PMMA (e.g. Telio® CAD A16)

- In preparation for adhesive cementation, apply an MMA-based light-curing conditioner on the bonding surface of the Telio CAD structure.

For more detailed information please also refer to the Instructions for Use of the conditioner.

2.3 Cementation with Multilink Hybrid Abutment

- Have the cleaned and conditioned components (ceramic structure, titanium bonding base) ready at hand for cementation.
- **The subsequent cementation procedure must be carried out quickly and without interruption. The working time of Multilink Hybrid Abutment is approx. 2 min at 23 °C / 73 °F (± 1 °C / 2 °F).**
- As a general rule, a new mixing tip is attached to the Multilink Hybrid Abutment syringe prior to each use.
- Apply a thin layer of Multilink Hybrid Abutment directly from the mixing tip to the bonding surface of the bonding base and **the bonding surface of the ceramic or PMMA structure.**
- Leave the mixing tip on the automix syringe until the next use. The remaining cement polymerizes in the tip and functions as a seal.
- Place the structure on the base in such a way that the position markings are aligned.
- Press the parts lightly and evenly together and check the correct relative position of the components (transition bonding base/structure).
- Subsequently, tightly press the components together for 5 s.
- Carefully remove excess in the screw channel, e.g. with a microbrush or brush, using rotary movements.
- **Important: Do not remove circular excess cement before curing has started, i.e. 3 minutes after mixing. Use a suitable instrument for this purpose (e.g. Le Cron). Secure the components in place using light pressure.**
- Apply glycerine gel (e.g. Liquid Strip) to the cementation joint to prevent the formation of an inhibition layer. The glycerine gel must be applied cautiously to avoid blending it with or displacing the composite. The gel must be left on the cementation joint until polymerization is complete.
- Next, the composite cement is left to auto-polymerize completely within 7 min.
- **Important: Do not move the components until Multilink Hybrid Abutment has completely cured. Secure them in place using e.g. diamond-coated tweezers.**
- After completion of the auto-polymerization, rinse off the glycerine gel with water.
- **Cautiously polish the cementation joint with rubber polishers at low speed (< 5,000 rpm) to avoid overheating**
- Remove any remaining cement residue in the screw channel with suitable rotating instruments.
- Steam clean the restoration.

3 Safety notes

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclar.com and your responsible competent authority.
- The current Instructions for Use are available in the download section of the Ivoclar Vivadent AG website (www.ivoclar.com).
- Explanation of the symbols: www.ivoclar.com/elfu
- The Summary of Safety and Clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>.
- Basic UDI-DI: 76152082ACEME003DZ

Warnings

- Observe the Safety Data Sheet (SDS) (available in the download section of the Ivoclar Vivadent AG website www.ivoclar.com).
- Prevent any contact of uncured materials with the skin/mucous membrane and eyes.
- Uncured material may cause slight irritation and may lead to a sensitization against methacrylates.
- Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Disposal information

Remaining stocks must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks.

The following clinical residual risks are known:

- Failure of the adhesive bond

4 Shelf life and storage

- Storage temperature 2–28 °C (36–82 °F)
- Store automix syringe with mixing tip attached after use.
- Do not use the product after the indicated date of expiration.
- Expiry date: see information on syringes and packages.

Before use, visually inspect the packaging and the product for damage. In case of any doubts, please contact Ivoclar Vivadent AG or your local dealer.

5 Additional information

Keep material out of the reach of children!

Not all products are available in all countries.

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions.

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