Tetric® N-Ceram 2

[en] Instructions for Use

Polymer-based dental restorative material
(intraoral light-curing)

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Patient target group

- Patients with permanent teeth
- Patients with deciduous teeth

Intended users / Special training

- Dentists
- No special training required.

Use

For dental use only.

Description

Tetric® N-Ceram 2 is a light-curing, radiopaque composite (300% Al) for the direct restorative treatment of anterior and posterior teeth (according to ISO 4049:2019 Type 1, Class 2, Group 1). Tetric N-Ceram 2 is also suitable for restoring occlusal surfaces.

100% aluminium has a radiopacity equivalent to that of dentin and 200% aluminium is equivalent to enamel.

Tetric N-Ceram 2 cures with light in the wavelength range of 400–500 nm.

Indications

Missing tooth structure in anterior and posterior teeth

Types of restorations

NOTE: Please observe the limitations of use.

- Restoration of permanent teeth
- Masking of discoloured tooth structure in permanent anterior teeth
- Splinting of mobile permanent teeth
- Repair of composite/ceramic restorations in permanent teeth
- Restoration of deciduous teeth. The limitations of use must be observed.

Contraindications

The use of the product is contraindicated if the patient is known to be allergic to any of its ingredients.

Limitations of use

- If a dry working field cannot be established.
- If the stipulated working procedures cannot be applied.
- Light intensities $>1,300 \text{ mW/cm}^2$ must not be used for deciduous teeth.
-  If the material is applied directly into the oral cavity of the patient, the Cavifil must not be used for more than one patient due to hygienic reasons (prevention of cross-contamination between patients).

Side effects

In individual cases, components of Tetric N-Ceram 2 may lead to sensitization. The product must not be used in such cases. In order to avoid irritation of the pulp, provide areas close to the pulp with suitable pulp/dentin protection. Selectively apply a calcium hydroxide-containing material in areas close to the pulp and cover with a suitable cavity liner.

Interactions

Substances containing eugenol/clove oil may inhibit the polymerization of methacrylate-based materials. Consequently, application of such materials together with Tetric N-Ceram 2 must be avoided. Discolouration may occur in combination with cationic mouthwashes, plaque disclosing agents and chlorhexidine.

Clinical benefit

- Reconstruction of chewing function
- Restoration of esthetics

Composition

Barium glass, copolymer, UDMA, Si-Zr mixed oxide, ytterbium trifluoride, Bis-EMA, Bis-GMA

Total content of inorganic fillers: 54–56 vol%

Particle size of inorganic fillers: between $0.11 \mu\text{m}$ and $15.5 \mu\text{m}$.

2 Usage

I. Shade determination

Clean the teeth prior to shade determination. The shade is selected with the tooth still moist using a shade guide (e.g. the Tetric N-Ceram 2 shade guide).

II. Isolation

Adequate relative or absolute isolation is required.

III. Cavity preparation

The cavity is prepared according to the principles of the adhesive technique, i.e. by preserving as much of the tooth structure as possible. Do not prepare sharp, internal edges and angles. Do not prepare additional undercuts in caries-free areas. The dimensions of the cavity are generally determined by the extent of the caries or the size of the old restoration. Bevel the enamel margins in the anterior region. Lightly round any sharp enamel edges in the posterior region (finishing diamonds, $25\text{--}40 \mu\text{m}$). Caries-free cervical defects are not prepared, only cleaned with pumice or other suitable cleaning pastes with the help of rubber cups or rotary brushes. Remove all residue in the cavity with water spray. Dry the cavity with water- and oil-free air.

IV. Pulp protection / Base

Do not apply a base material when using an enamel/dentin bonding agent. In very deep cavities only, areas close to the pulp should be selectively coated with a calcium hydroxide liner. Subsequently cover with a pressure resistant cement. Do not cover the remaining cavity walls so that they can be used to generate a bond with an enamel-dentin adhesive.

V. Placement of matrix / interdental wedge

Use a wrap-around matrix for cavities affecting the proximal area or a sectional matrix band and secure it with wedges.

VI. Conditioning / Application of the bonding agent

Apply the bonding agent to the tooth structure according to the instructions for use of the manufacturer of the product used. Ivoclar recommends using a universal adhesive.

VII. Application of Tetric N-Ceram 2

- For optimum results, apply Tetric N-Ceram 2 in layers of max. 2 mm thickness or 1.5 mm thickness (Dentin shades) using Cavifil Injector or a comparable applicator and adapt it with a suitable instrument.
- Adapt the material correctly in order to ensure an intimate contact of the composite resin with the cavity walls.
- Prevent incomplete polymerization of the restoration by ensuring sufficient exposure to the curing light.
- For the recommendations regarding exposure time (Exposure time) per increment and light intensity (Light intensity) see table 1 (Table 1).
-  The instructions for use of the curing light must be observed.

- When using a metal matrix, additionally polymerize the composite material from the buccal and lingual/palatal aspect after removing the matrix, if no Bluephase® N curing light is used.
- If the light guide could not be positioned ideally, e.g. at a distance to the composite or at a divergent light emission angle, light-cure the composite material again.
- If the diameter of the light guide tip is smaller than the restoration diameter, perform an overlapping polymerization in order to ensure that all areas of the restoration are covered.
- Optionally, a flowable composite can be applied as an initial layer.
 Cure this layer separately according to the respective instructions for use.
- Application in the case of repairs: Roughen the surface to be repaired using diamond grinders and thoroughly clean with water spray and dry with oil- and water-free compressed air. The surface of ceramic restorations must be treated with a suitable ceramic primer according to the manufacturer's instructions.
 Apply the bonding agent to the composite / ceramic surface according to the instructions for use of the product in use. Apply the composite resin as described in section VII.

VIII. Finishing / Checking the occlusion / Polishing

Remove any excess material with tungsten carbide or diamond finishers after polymerization. Check the occlusion and articulation and apply appropriate grinding corrections to prevent premature contacts or undesired articulation paths on the surface of the restoration. Use polishers as well as polishing discs and polishing strips to polish the restoration to a high gloss.

Application notes

- In the case of repairs, additional Tetric N-Ceram 2 can be directly applied to polymerized material. If the Tetric N-Ceram 2 restoration has already been polished, it must first be roughened and wetted with an adhesive before fresh material is applied.
- Tetric N-Ceram 2 should have ambient temperature when applied. Cool temperatures render the material difficult to extrude.
- Syringes or Cavifils should not be disinfected with oxidizing disinfection agents.

3 Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclar.com, and your responsible competent authority.
- The current Instructions for Use are available on the website: www.ivoclar.com.
- Explanation of symbols: www.ivoclar.com/elFU
- The Summary of Safety and Clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>.
- Basic UDI-DI: 76152082AFILL004JG

Warnings

- Observe the Safety Data Sheet (SDS) (available at www.ivoclar.com).
- Unpolymerized Tetric N-Ceram 2 should not come in contact with skin, mucous membrane, or eyes. Unpolymerized Tetric N-Ceram 2 may have a slight irritating effect and may lead to a sensitization against methacrylates. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Disposal information

Remaining stocks or removed restorations must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks.

The following clinical residual risks are known:

- Failure of the adhesive bond (loss of the filling)
- Postoperative sensitivity
- Inclusion of air bubbles during the placement of the filling
- Chipping, fractures
- Ingestion of material

4 Shelf life and storage

- Storage temperature: 2 - 28 °C (36 - 82 °F)
- Close syringes/Cavifils immediately after usage. Exposure to light leads to premature polymerization.
- Do not use the product after the indicated date of expiration.
- Expiry date: see information on Cavifils, syringes and packages.
- Before use, visually inspect the packaging and the product for damage. In case of any doubts, please contact Ivoclar Vivadent AG or your local dealer.

5 Additional information

Keep material out of the reach of children!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the product for its suitability and use for any purpose not explicitly stated in the Instructions.

Table 1

Light intensity	Exposure time
500 – 900 mW/cm ²	20 s
1,000 – 1,300 mW/cm ²	10 s
1,800 – 2,200 mW/cm ²	5 s

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