Tetric® N-PowerFill 2

[en] Instructions for Use

Polymer-based dental restorative material
(intraoral light-curing)

Rx ONLY

C € 0123

Date information prepared:
2023-09-13/Rev. 1

Manufacturer:
Ivoclar Vivadent AG
Bendererstrasse 2
9494 Schaan / Liechtenstein
www.ivoclar.com

ivoclar

English

1 Intended use

Intended purpose

Direct restorations of posterior teeth

Patient target group

- Patients with permanent teeth
- Patients with deciduous teeth

Intended users / Special training

- Dentists
- No special training required.

Use

For dental use only.

Description

Tetric® N-PowerFill 2 is a light-curing, radiopaque composite (300% AI) for the direct restorative treatment of posterior teeth (according to ISO 4049:2019 Type 1, Class 2, Group 1). Tetric N-PowerFill 2 is also suitable for restoring occlusal surfaces. 100% aluminium has a radiopacity equivalent to that of dentin and 200% aluminium is equivalent to enamel.

Tetric N-PowerFill 2 cures with light in the wavelength range of 400–500 nm and can be applied in layers of up to 4 mm.

Indications

Missing tooth structure in posterior teeth (Classes I and II)

Types of restorations:

- Restorations in permanent posterior teeth (Classes I and II)
- Restorations in posterior deciduous teeth (Class I and II).
The limitations of use must be observed.
- Reconstructive build-ups, when light-cured at a light intensity of $\leq 1,300 \text{ mW/cm}^2$.

Contraindications

- The use of the product is contraindicated if the patient is known to be allergic to any of its ingredients.

Limitations of use

- Due to esthetic reasons, this product is not suitable for Class III and IV restorations.
- Avoid direct exposure of the gingiva, mucous membrane or skin when using the 3sCure mode.
- The 3sCure mode must not be used in case of caries profunda and very deep cavities.
- If a dry working field cannot be established.
- If the stipulated working procedures cannot be applied.

- The light intensity $>1300 \text{ mW/cm}^2$ must not be used for deciduous teeth.
- When light-cured with the Bluephase PowerCure in the 3sCure curing mode, Tetric N-PowerFill 2 can be used for: restorations in the posterior region of permanent dentition (Classes I and II) when light-cured from the occlusal aspect.
- Apply the product at ambient temperature. Cool temperatures render the material difficult to extrude.

Side effects

In rare cases, components of Tetric N-PowerFill 2 may lead to sensitization. The product must not be used in such cases. In order to avoid irritation of the pulp, provide areas close to the pulp with suitable pulp / dentin protection. Selectively apply a calcium hydroxide-containing material in areas close to the pulp and cover with a suitable cavity liner.

Interactions

Phenolic substances such as eugenol / clove oil inhibit the polymerization of methacrylate-based materials. Consequently, the application of such materials together with Tetric N-PowerFill 2 must be avoided. Discolouration may occur in combination with cationic mouthwashes, plaque disclosing agents and chlorhexidine. Disinfectants with an oxidative effect (e.g. hydrogen peroxide) may interact with the initiator system, which in turn may impair the curing process.

Clinical benefit

- Reconstruction of chewing function
- Restoration of esthetics

Composition

Barium glass, copolymer, Si-Zr mixed oxide, Bis-GMA, ytterbium trifluoride, Bis-PMA, UDMA, Bis-EMA

Total content of inorganic fillers: 53–54 vol%

Particle size of the inorganic fillers: between $0.11 \mu\text{m}$ and $15.46 \mu\text{m}$

2 Application

I. Shade determination

Clean the teeth prior to shade determination. The shade is selected with the tooth still moist using a shade guide (e.g. the Tetric N-PowerFill/Flow 2 shade guide).

II. Isolation

Adequate relative or absolute isolation is required.

III. Cavity preparation

The cavity is prepared according to the principles of the adhesive technique, i.e. by preserving as much of the tooth structure as possible. Do not prepare sharp, internal edges and angles. Do not prepare additional undercuts in caries-free areas. The dimensions of the cavity are generally determined by the extent of the caries or the size of the old restoration. Lightly round any sharp enamel edges in the posterior region (finishing diamonds, 25–40 μm). Caries-free cervical defects are not prepared, only cleaned with pumice or other suitable cleaning pastes with the help of rubber cups or rotary brushes. Remove all residue in the cavity with water spray. Dry the cavity with water- and oil-free air.

IV. Pulp protection / Base

Do not apply a base material when using an enamel-dentin bonding agent. In very deep cavities only, areas close to the pulp should be selectively coated with a calcium hydroxide liner. Subsequently cover with a pressure resistant cement. Do not cover the remaining cavity walls so that they can be used to generate a bond with an enamel-dentin adhesive.

V. Placement of matrix / interdental wedge

Use a wrap-around matrix or a sectional matrix band for cavities affecting the proximal area and secure it with wedges.

VI. Conditioning / Application of the bonding agent

Apply the bonding agent to the tooth structure according to the instructions for use of the manufacturer of the product used. Ivoclar recommends using a universal adhesive.

VII. Application of Tetric N-PowerFill 2

- In order to achieve optimum results, Tetric N-PowerFill 2 should be applied in increments of max. 4 mm and adapted to the cavity walls with a suitable instrument.

- Adapt the material correctly in order to ensure an intimate contact of the composite resin with the cavity walls.
- Prevent incomplete polymerization of the restoration by ensuring sufficient exposure to the curing light.
- For the recommendations regarding exposure time (Exposure time) per increment and light intensity (Light intensity) see table 1 (Table 1).
-  **The instructions for use of the curing light must be observed.**
- When using a metal matrix, additionally polymerize the composite material from the buccal and lingual/palatal aspect after removing the matrix, if no Bluephase® N curing light is used.
- If the light guide could not be positioned ideally, e.g. at a distance to the composite or at a divergent light emission angle, light-cure the composite material again.
- If the diameter of the light guide tip is smaller than the restoration diameter, perform an overlapping polymerization in order to ensure that all areas of the restoration are covered.
-  Optionally, a flowable composite can be applied as an initial layer. Cure this layer separately according to the respective instructions for use.

VIII. Finishing / Checking the occlusion / Polishing

Remove any excess material with tungsten carbide or diamond finishers after polymerization. Check the occlusion and articulation and apply appropriate grinding corrections to prevent premature contacts or undesired articulation paths on the surface of the restoration. Use polishers as well as polishing discs and polishing strips to polish the restoration to a high gloss.

Application notes

- Avoid direct exposure of the gingiva, mucous membrane or skin to the light emitted by the curing light.
- In the case of repairs, additional Tetric N-PowerFill 2 can be directly applied to the polymerized material. If the Tetric N-PowerFill 2 restoration has already been polished, it must first be roughened and wetted with an adhesive before fresh Tetric N-PowerFill 2 is applied.
- The recommended increment thickness is based on hardness profile measurements.
- Syringes must not be disinfected with oxidizing disinfection agents.
- Dispose of contaminated syringes.

3 Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclar.com, and your responsible competent authority.
- The current Instructions for Use are available on the website (www.ivoclar.com).
- Explanation of the symbols: www.ivoclar.com/elFU
- The Summary of Safety and Clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>.
- Basic UDI-DI: 76152082AFILL001JA

Warnings

- Observe the Safety Data Sheet (SDS) (available at www.ivoclar.com).
- Unpolymerized Tetric N-PowerFill 2 should not come in contact with skin, mucous membrane and eyes. Unpolymerized Tetric N-PowerFill 2 may have a slight irritating effect and may lead to a sensitization against methacrylates. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Disposal information

Remaining stocks or removed restorations must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks.

The following clinical residual risks are known:

- Failure of the adhesive bond (loss of the filling)
- Postoperative sensitivity
- Inclusion of air bubbles during the placement of the filling
- Wear of the filling
- Heat development during the curing procedure
- Chipping, fractures
- Ingestion of material

4 Shelf life and storage

- Storage temperature: 2–28 °C (36–82 °F)
- Close syringes immediately after usage. Exposure to light leads to premature polymerization.
- Do not use the product after the indicated expiry date.
- Expiry date: see information on syringes and packages.
- Before use, visually inspect the packaging and the product for damage. In case of any doubts, please contact Ivoclar Vivadent AG or your local dealer.

5 Additional information

Keep material out of the reach of children!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

Table 1

Light intensity	Exposure time
500 – 900 mW/cm ²	20 s
1,000 – 1,300 mW/cm ²	10 s
1,800 – 2,200 mW/cm ²	5 s
2,700 – 3,300 mW/cm ²	3 s

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