IPS Empress® Direct

[en] Instructions for Use

Polymer-based dental restorative material
(intraoral light-curing)

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For dental use only.

Description

IPS Empress Direct is a light-curing, radiopaque composite (210% AI) for the direct restorative treatment of anterior and posterior teeth (according to ISO 4049:2019 Type 1, Class 2, Group 1).

IPS Empress Direct is also suitable for restoring occlusal surfaces. 100% aluminum has a radiopacity equivalent to that of dentin and 200% aluminum is equivalent to enamel.

IPS Empress Direct cures with light in the wavelength range of 400–500 nm.

IPS Empress Direct Opal shade is classified as Type 2, Class 2, Group 1 in accordance with ISO 4049:2019. The IPS Empress Direct Opal shade is not suitable for the reconstruction of occlusal surfaces and it is not radiopaque.

Indications

Missing tooth structure in anterior and posterior teeth

Types of restorations:

- Fillings in permanent teeth
- Corrections of tooth positions and tooth shape
- Direct veneers
- Repair of composite restorations

Contraindications

The use of these products is contraindicated if the patient is known to be allergic to any of its ingredients.

Limitations of use

- If a dry working field cannot be established.
- If the stipulated working technique cannot be applied.
- Apply the product at ambient temperature. Cool temperatures render the material difficult to extrude.
-  If IPS Empress Direct is applied directly in the oral cavity of the patient, the Cavifil must not be used for more than one patient due to hygienic reasons (prevention of cross-contamination between patients).

Side effects

In rare cases, components of IPS Empress Direct may lead to sensitization. The products must not be used in such cases.

In order to avoid irritation of the pulp, provide areas close to the pulp with suitable pulp/dentin protection. Selectively apply a calcium hydroxide-containing material in areas close to the pulp and cover with a suitable cavity liner.

Interactions

- Phenolic substances such as eugenol/clove oil inhibit the polymerization of methacrylate-based materials. Consequently, the application of such materials together with IPS Empress Direct must be avoided.
- Discoloration may occur in combination with cationic mouthwashes, plaque disclosing agents and chlorhexidine.
- Disinfectants with an oxidative effect (e.g. hydrogen peroxide) may interact with the initiator system, which in turn may impair the curing process.

Clinical benefit

- Restoration of esthetics
- Reconstruction of chewing function

Composition

- **IPS Empress Direct, Enamel Shades:**
Barium glass, Si-Zr mixed oxide, UDMA, DCP, Bis-PMA, Bis-GMA, barium aluminum fluorosilicate glass
- **IPS Empress Direct, Dentin Shades:**
Barium glass, co-polymer, UDMA, Si-Zr mixed oxide, ytterbium trifluoride, DCP, Bis-PMA, Bis-GMA, barium aluminum fluorosilicate glass
- **IPS Empress Direct, Trans 20/Trans 30:**
barium glass, UDMA, Si-Zr mixed oxide, co-polymer, DCP, Bis-PMA, Bis-GMA
- **IPS Empress Direct, Opal:**
Co-polymer, silicone oxide, aromatic-aliphatic UDMA, D3MA, DCP, UDMA

Total content of inorganic fillers: 52–59 vol.% (opalescent shade: 45 vol.%).
Particle size of inorganic fillers: between 0.03 µm and 16.3 µm (opalescent shade: 0.14 µm–13.8 µm).

2 Application

Shade determination

Clean the teeth prior to shade determination. The shade is selected with the tooth still moist using a shade guide (e.g. the IPS Empress Direct shade guide).

Isolation

Adequate relative or absolute isolation is required.

Cavity preparation

The cavity is prepared according to the principles of the adhesive technique, i.e. by preserving as much of the tooth structure as possible. Do not prepare sharp internal edges and angles. Do not prepare additional undercuts in caries-free areas. The dimensions of the cavity are generally determined by the extent of the caries or the size of the old restoration. Bevel the enamel margins in the anterior region. Lightly round any sharp enamel edges in the posterior region (finishing diamonds, 25–40 µm). Caries-free cervical defects are not prepared, but only cleaned with pumice or other suitable cleaning pastes with the help of rubber cups or rotary brushes. Remove all residue in the cavity with water spray. Dry the cavity with water- and oil-free air.

Pulp protection/Base

Do not apply a base material when using an enamel-dentin bonding agent. In very deep cavities, areas close to the pulp should be selectively coated with a calcium hydroxide liner. Then, cover with a pressure-resistant cement. Do not cover the remaining cavity walls so that they can be used to generate a bond with an enamel-dentin adhesive.

Placement of matrix/interdental wedge

Use a wrap-around matrix or a sectional matrix band for cavities affecting the proximal area and secure it with wedges.

Conditioning/Application of the bonding agent

Apply the bonding agent to the tooth structure according to the instructions for use of the manufacturer of the product used. Ivoclar recommends using a universal adhesive.

Application of the product:

- For optimum results, apply IPS Empress Direct in layers of max. 2 mm thickness or 1.5 mm thickness (Dentin shade) using a suitable instrument. When Cavifils are used, apply the material from the Cavifil with the Cavifil injector or with a comparable applicator.
- Adapt the material so that an intimate contact of the composite resin with the cavity walls is ensured.
- Prevent incomplete polymerization of the restoration by ensuring sufficient exposure to the curing light.
- For the recommendations regarding exposure time per increment (increment thickness) and light intensity see Table 1.

NOTICE! The instructions for use of the curing light used must be observed.

- When using a metal matrix, additionally polymerize the composite material from the buccal and lingual/palatal aspect after removing the matrix, if no Bluephase® curing light is used.
- If the light guide could not be positioned ideally, e.g. at a distance to the composite or at a divergent light emission angle, light-cure the composite material again.
- If the diameter of the light guide tip is smaller than the restoration diameter, perform an overlapping polymerization in order to ensure that all areas of the restoration are covered.

Observe the instructions for use of the product used.

- Optionally, a flowable composite can be applied as an initial layer. Cure this layer separately.

Observe the instructions for use of the product used.

- Application in the case of repairs:
Roughen the surface to be repaired using diamond burs and thoroughly clean with water spray and dry with oil- and water-free compressed air.
Apply the bonding agent to the composite surface according to the instructions for use of the product used. Apply the composite resin as described in section VII.

Finishing/Checking the occlusion/Polishing

After polymerization, remove any excess material with tungsten carbide or diamond finishers. Check the occlusion and articulation and apply appropriate grinding corrections to prevent premature contacts or undesired articulation paths on the surface of the restoration. Use polishers as well as polishing discs and polishing strips to polish the restoration to a high gloss.

Application notes

- If necessary, the products can also be applied directly to already polymerized material. If the restoration is already polished, it must be roughened and wetted first with an adhesive before the fresh product can be applied.
(IPS Empress Direct Opal shade is not suitable for the reconstruction of occlusal surfaces).
- Do not expose the product to intensive light while applying it, as this might lead to premature polymerization.
- Syringes or Cavifils must not be disinfected with oxidizing disinfection agents.
- Avoid any contact of the reusable syringe not covered with a protective sleeve with the mouth of the patient.
- Dispose of contaminated syringes.

3 Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclar.com, and your responsible competent authority.
- The current Instructions for Use are available on the website (www.ivoclar.com/eIFU).
- Explanation of the symbols: www.ivoclar.com/eIFU
- The Summary of Safety and Clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>.
- Basic UDI-DI: 76152082AFILL009JS

Warnings

- Observe the Safety Data Sheet (SDS) (available at www.ivoclar.com).
- Avoid contact of the uncured product with the skin/mucous membrane or eyes. The unpolymerized product may have a slight irritating effect and may lead to a sensitization against methacrylates. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Disposal information

Remaining stocks or removed restorations must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks.

The following clinical residual risks are known:

- Failure of the adhesive bond (loss of the filling)
- Postoperative sensitivity
- Inclusion of air bubbles during the placement of the filling
- Chipping, fractures
- Ingestion of material

4 Shelf life and storage

- Storage temperature 2–28 °C (36–82 °F)
- Close syringes/Cavifils immediately after usage. Exposure to light leads to premature polymerization.
- Do not use the product after the indicated date of expiration.
- Expiry date: See information on Cavifils, syringes and packages.
- Before use, visually inspect the packaging and the product for damage. In case of any doubts, please contact Ivoclar Vivadent AG or your local dealer.

5 Additional information

Keep material out of the reach of children!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the materials for their suitability and use for any purpose not explicitly stated in the Instructions.

Table 1

Light intensity	Exposure time
500–900 mW/cm ²	20 s
1,000–1,300 mW/cm ²	10 s

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