Tetric® N-Effect 2

[en] Instructions for Use

Polymer-based dental restorative material
(intraoral light-curing)

C € 0123

Manufacturer:
Ivoclar Vivadent AG
Bendererstrasse 2
9494 Schaan/Liechtenstein
www.ivoclar.com

Rx ONLY

Date information prepared:
2024-08-12 / Rev. 0
web

ivoclar

English

1 Intended use

Intended purpose

Esthetic characterization of direct and indirect anterior and posterior restorations and masking of discolorations.

Patient target group

Patients with permanent teeth.

Intended users

Dentists

Special training

No special training required.

Use

For dental use only.

Description

Tetric N-Effect 2 is a flowable, light-curing composite for direct restorative procedures, which is used to characterize esthetically demanding composite restorations (not suitable for the reconstruction of occlusal surfaces – according to ISO 4049:2019 Type 2, Class 2, Group 1).

Tetric N-Effect 2 cures with light in the wavelength range of 400–500 nm.

Indications

None

Types of restorations:

- As an intermediate or covering layer in the fabrication of esthetically demanding composite restorations.
- Repair of composite restorations.

Contraindications

The use of these products is contraindicated if the patient is known to be allergic to any of its ingredients.

Limitations of use

- If a dry working field cannot be established.
- If the stipulated working technique cannot be applied.
- Apply the product at ambient temperature. Cool temperatures render the material difficult to extrude.
- The use of different application cannulas can render the material difficult to extrude.
- The syringe may be used for multiple applications, but should be protected by a single-use hygienic sleeve. Use a suitable hygienic protective sleeve for the syringe as indicated by the manufacturer.
-  The application cannula of the syringe must not be used for more than one patient due to hygienic reasons (prevention of cross-contamination between patients).

Side effects

In rare cases, components of Tetric N-Effect 2 may lead to sensitization. The products must not be used in such cases.

In order to avoid irritation of the pulp, provide areas close to the pulp with suitable pulp/dentin protection. Selectively apply a calcium hydroxide-containing material in areas close to the pulp and cover with a suitable cavity liner.

Interactions

- Phenolic substances such as eugenol/clove oil inhibit the polymerization of methacrylate-based materials. Consequently, the application of such materials together with Tetric N-Effect 2 must be avoided.
- Discoloration may occur in combination with cationic mouthwashes, plaque disclosing agents and chlorhexidine.
- Disinfectants with an oxidative effect (e.g. hydrogen peroxide) may interact with the initiator system, which in turn may impair the curing process.

Clinical benefit

Restoration of esthetics

Composition

Co-polymer, aromatic-aliphatic UDMA, silanized silicone oxide, D3MA, DCP, UDMA

Total content of inorganic fillers: 34 vol%.

Particle size of inorganic fillers: between 0.14 µm and 13.8 µm.

2 Application

Shade determination

The product is available in the shade trans opal.

Isolation

Adequate relative or absolute isolation is required.

Cavity preparation

The cavity is prepared according to the principles of the adhesive technique, i.e. by preserving as much of the tooth structure as possible. Do not prepare sharp internal edges and angles. Do not prepare additional undercuts in caries-free areas. The dimensions of the cavity are generally determined by the extent of the caries or the size of the old restoration. Bevel the enamel margins in the anterior region. Lightly round any sharp enamel edges in the posterior region (finishing diamonds, 25–40 µm). Caries-free cervical defects are not prepared, but only cleaned with pumice or other suitable cleaning pastes with the help of rubber cups or rotary brushes. Remove all residue in the cavity with water spray. Dry the cavity with water- and oil-free air.

Placement of matrix/interdental wedge

Use a wrap-around matrix or a sectional matrix band for cavities affecting the proximal area and secure it with wedges.

Conditioning/Application of the bonding agent

Apply the bonding agent to the tooth structure according to the instructions for use of the manufacturer of the product used. Ivoclar recommends using a universal adhesive.

Application of the product:

- For optimum results, apply Tetric N-Effect 2 in layers of max. 2 mm thickness.
 - Keep the tip of the application cannula immersed in the material during the application of Tetric N-Effect 2 to avoid the formation of air entrapments.
 - Adapt the material so that an intimate contact of the composite resin with the cavity walls is ensured, using a suitable instrument (e.g. a probe).
 - Prevent incomplete polymerization of the restoration by ensuring sufficient exposure to the curing light.
 - For the recommendations regarding exposure time per increment (increment thickness) and light intensity see Table 1.
-  **NOTICE! The instructions for use of the curing light used must be observed.**
- When using a metal matrix, additionally polymerize the composite material from the buccal and lingual/palatal aspect after removing the matrix, if no Bluephase® N curing light is used.
 - If the light guide could not be positioned ideally, e.g. at a distance to the composite or at a divergent light emission angle, light-cure the composite material again.

- If the diameter of the light guide tip is smaller than the restoration diameter, perform an overlapping polymerization in order to ensure that all areas of the restoration are covered.
 Observe the instructions for use of the product used.
- Optionally, a flowable composite can be applied as an initial layer. Cure this layer separately.
 Observe the instructions for use of the product used.
- Application in the case of repairs:
Roughen the surface to be repaired using diamond burs and thoroughly clean with water spray and dry with oil- and water-free compressed air.
Apply the bonding agent to the composite surface according to the instructions for use of the product used. Apply the composite resin as described in section VII.

Finishing/Checking the occlusion/Polishing

After polymerization, remove any excess material with tungsten carbide or diamond finishers. Check the occlusion and articulation and apply appropriate grinding corrections to prevent premature contacts or undesired articulation paths on the surface of the restoration. Use polishers as well as polishing discs and polishing strips to polish the restoration to a high gloss.

Application notes

- If necessary, the product can also be applied directly to already polymerized material. If the restoration is already polished, it must be roughened and wetted first with an adhesive before the fresh product can be applied.
(Tetric N-Effect 2 is not suitable for the reconstruction of occlusal surfaces).
- Do not expose the product to intensive light while applying it, as this might lead to premature polymerization.
- Syringes must not be disinfected with oxidizing disinfection agents.
- Use a suitable hygienic protective sleeve for the syringe as indicated by the manufacturer.
- Avoid any contact of the reusable syringe not covered with a protective sleeve with the mouth of the patient.
- Dispose of contaminated syringes.

3 Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclar.com, and your responsible competent authority.
- The current Instructions for Use are available on the website (www.ivoclar.com/eIFU).
- Explanation of the symbols: www.ivoclar.com/eIFU
- The Summary of Safety and Clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>.
- Basic UDI-DI: 76152082AFILL015JM

Warnings

- Observe the Safety Data Sheet (SDS) (available at www.ivoclar.com).
- Avoid contact of the uncured product with the skin/mucous membrane or eyes. The unpolymerized product may have a slight irritating effect and may lead to a sensitization against methacrylates. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Disposal information

Remaining stocks or removed restorations must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks.

The following clinical residual risks are known:

- Failure of the adhesive bond (loss of the filling)
- Postoperative sensitivity
- Inclusion of air bubbles during the placement of the filling
- Chipping, fractures
- Ingestion of material

4 Shelf life and storage

- Storage temperature 2–28 °C (36–82 °F)
- Close syringes immediately after usage. Exposure to light leads to premature polymerization.
- Do not use the product after the indicated date of expiration.
- Expiry date: See information on syringes and packages.
- Before use, visually inspect the packaging and the product for damage. In case of any doubts, please contact Ivoclar Vivadent AG or your local dealer.

5 Additional information

Keep material out of the reach of children!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the materials for their suitability and use for any purpose not explicitly stated in the Instructions.

Table 1

Light intensity	Exposure time
500–900 mW/cm ²	20 s
1,000–1,300 mW/cm ²	10 s

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