Intended use

Intended purpose

Absolute isolation

Use

For dental use only.

Description

OptraDam® Plus is an anatomically shaped rubber dam for the absolute isolation of the working field. It provides retraction of lips and cheeks during dental treatment and offers the clinician increased visibility and accessibility. In addition, the device facilitates the evacuation of water spray. OptraDam Plus features high flexibility and elasticity in all directions of movement, provides patient comfort and assists in keeping the patient's mouth open. OptraDam Plus is available in regular and small sizes.

It is suitable for use in therapeutic dental procedures, for which an absolutely dry and isolated working field is required (e.g. for the adhesive cementation of restorations, direct restorative procedures and root canal treatments).

Indication

Isolation of single and multiple teeth during dental treatment

Contraindications

The use of OptraDam Plus is contraindicated in patients with latex allergy.

Limitations of use

⊗ OptraDam Plus is for single patient use only. The product is not intended to be reprocessed or reused. The material is not suitable for sterilization. Cleaning and disinfecting agents may have a detrimental effect on the material and lead to undesirable reactions in the patient.

Side effects

 OptraDam Plus contains natural rubber latex which may cause allergic reactions.

Composition

OptraDam Plus is made of latex and POM and it is powdered with corn starch.

Application

When OptraDam Plus is used, the small intraoral ring of the device is positioned in the area of the gingivobuccal fold, while the outer ring remains outside the mouth.

The elastic component between the two rings thus embraces the lips of the patient and provides retraction due to the restoring force of the rings (see Figure 1).

# OptraDam® Plus

EN Instructions for Use – Anatomically shaped rubber dam	NO Bruksanvisning – Anatomisk formet kofferdam	LV Lietošanas instrukcija – Anatomiskas formas gumijas koferdams
DE Gebrauchsinformation – Anatomisch ausgeformter Kofferdam	NL Gebruiksaanwijzing – Anatomisch gevormde cofferdam	LT Naudojimo informacija – Anatomiška is uformuota guminė užtvara
FR Mode d'emploi – Digue anatomique en latex	EL Οδηγίες Χρήσεως – Ανατομικά διαμορφωμένος ελαστικός απομονωτήρας	SK Návod na používanie – Anatomicky tvarovaný koferdam
IT Istruzioni d'uso – Diga di gomma anatomica	TR Kullanma Talimatı – Anatomik diarak şekillendirilmiş rubber dam	HU Használati utasítás – Anatómiailag formált koferdám
ES Instrucciones de uso – Dique de goma anatómico	RU Инструкция по применению – Коффердам анатомической формы	SR Упутство за употребу – Анатомски обликовангуменикофердам
PT Instruções de Uso – Dique de borracha com conformação anatómica	PL Instrukcja stosowania – Anatomycznie ukształtowany koferdam	MK Упатство за употреба – Анатомски обликуван кофердам
SV Bruksanvisning – Anatomiskt utformat kofferdam	SL Navodila za uporabo – AAнатомско oblikovana gumijasta opna	BG Инструкции за употреба – Анатомично оформен кофердам
		SQ Udhëzime përdorimi – Dige gome me formë anatomike
		RO Instrucțiuni de utilizare – Diga din cauciuc cu forma anatomica

Fig. 1: OpraDam Plus in place

1. Select the appropriate sized OpraDam Plus. The size (R for regular or S for small) is indicated on the packaging and on the back of the latex bag.
2. The rubber dam punch is used to make the appropriate holes, depending on the teeth to be exposed. Check the quality of your rubber dam punch regularly. Blunt punches perforate the latex only incompletely, which results in the latex material tearing and thus in insufficient isolation of the teeth.
3. OpraDam Plus may be inserted from both the 12 o'clock position with the patient lying down, or from the 8 o'clock position with the patient sitting upright (or, the 4 o'clock position if the clinician is left-handed).
4. To ease placement, the operator should grasp the intraoral ring by reaching through the outer ring. The intraoral ring should be held between the thumb and middle finger and pressed together slightly (see Figure 2). When grasping the intraoral ring, make sure that it is not completely pulled out of the latex bag.

Fig. 2: Correct hold on the intraoral ring with thumb and middle finger. The ring is slightly compressed in the process.

5. The slightly compressed intraoral ring is then inserted into the buccal corridor on one side, so that the elastic component embraces the corner of the mouth and the outer ring is still located extraorally (see Figure 2).
6. As soon as OpraDam Plus is securely in place on one side, the intraoral ring is inserted in the same manner on the other side by slightly bending it (see Figure 2).
7. Subsequently, the intraoral ring is placed behind the lower and upper lips and thus the OpraDam Plus assumes its final, stable position (see Figure 3). The intraoral ring can be inserted with more ease if the patient's mouth is relaxed. If, in individual cases, the intraoral ring should show a tendency to slip out of the gingivobuccal fold if the mouth is closed completely, positioning of the intraoral ring deeper inside the vestibule is usually sufficient to solve the problem. Sometimes, the use of another size may also be a solution.
8. The dam is secured around the teeth by fitting the septum of the dam interproximally and in the sulcus using dental floss.

Fig. 3: After having positioned the intraoral ring behind the corners of the mouth, it is pushed behind the lower and upper lip.

9. The securely positioned OpraDam Plus improves the overall view and accessibility of the working field. The saliva ejector is introduced below the integrated frame and positioned in the corner of the mouth. Lateral movements of the mandible and opening of the mouth to various degrees may provide additional space and access to the treatment area. If the interproximal spaces are large and retention of the rubber dam is insufficient, a ligature of dental floss may be used additionally to hold the latex bag in place.
10. Removal of the rubber dam: Carefully cut the interdental latex septa. Using a paper towel, grasp the outer ring in the lower area and lift it slightly so that the intraoral ring is loosened from the lower gingivobuccal fold. The intraoral ring can then be easily removed from the upper gingivobuccal fold and hygienically disposed of with the paper towel.

#### Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: [www.ivoclar.com](http://www.ivoclar.com), and your responsible competent authority.
- The current Instructions for Use are available in the download section of the Ivoclar Vivadent AG website ([www.ivoclar.com](http://www.ivoclar.com)).

#### Disposal information

- Remaining stocks must be disposed of according to the corresponding national legal requirements.

#### Shelf life and storage

- Storage temperature 2–28 °C (36–82 °F)
- Keep out of direct sunlight.
- Do not use the product after the indicated date of expiration.
- Expiration date: See note on packaging

#### Additional information

Keep out of the reach of children!

The material may be developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the instructions or the stipulated area of use. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

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