OptraGate®

[en] Instructions for Use

Flexible aid for full mouth access

Rx ONLY

Manufacturer:
Ivoclar Vivadent AG
Bendererstrasse 2
9494 Schaan/Liechtenstein
www.ivoclar.com

Date information prepared:
2024-11-27 / Rev. 2

ivoclar

English

1 Intended use

Intended purpose

Access to the oral cavity of the patient

Patient target group

- Patients with permanent teeth
- Patients with deciduous teeth
- Edentulous patients

Intended users

- Dentists (clinical procedure)
- Dental assistants (assistance during the clinical procedure)
- Dental hygienists (tooth cleaning procedure in the dental practice)

Special training

No further special training is needed.

Use

For dental use only.

Description

OptraGate® is an auxiliary aid for clinical use which facilitates access to the working field. It provides retraction of lips and cheeks during dental treatment, offers the operator increased visibility and accessibility and facilitates saliva and moisture control in the oral cavity. The device, which features three-dimensional flexibility and elasticity, ensures patient comfort and assists in keeping the patient's mouth open. It is available in sizes Regular, Small and Junior.

Indications

None

Areas of application:

The product is designed to retract the lips and cheeks during diagnostic, preventive, therapeutic and orthodontic dental procedures, e.g.

- initial examination;
- professional teeth cleaning;
- fissure sealing;
- tooth whitening;
- periodontal treatment;
- restorative procedures;
- preparation;
- impression taking;
- cementation of restorations and orthodontic retention elements.

Contraindications

The use of the product is contraindicated if the patient is known to be allergic to any of its ingredients.

Limitations of use

- OptraGate is not sterile and thus not suitable for use in sterile oral surgical or implantological procedures.
- OptraGate is not suitable for use in treatment procedures requiring absolute isolation.
-  For single use only.
- The product is not intended to be reprocessed or reused. The material is not suitable for sterilization. Cleaning and disinfecting agents may have a detrimental effect on the material and lead to undesirable reactions in the patient.

Side effects

There are no known side effects to date.

Interactions

There are no known interactions to date.

Clinical benefit

Easier access to the treatment field.

Composition

Thermoplastic elastomer (SEBS) and polypropylene. OptraGate is latex-free.

2 Application

The elastic material of OptraGate is stretched by means of two rings. The thicker ring of the device is positioned within the gingivobuccal fold, while the thinner ring, which features two tabs, remains outside the mouth. The elastic component between the two rings thus embraces the lips of the patient and provides retraction due to the restoring force of the rings (see Figure 1).

Fig 1: OptraGate in place with the tabs pointing downward. It does not interfere with the patient's occlusion.

1. Select the appropriate size of the OptraGate. The selection of the appropriate size should be guided by the distance between the corners of the mouth when the mouth is in a relaxed position (i.e. with the lips closed or slightly open). In order to obtain sufficient retraction of the oral soft tissues, the diameter of the inner ring should ideally be approx. 2 cm wider than the mouth corner distance.

Therefore the following guidelines apply:

- Size Junior: Mouth corner distance \approx 4,5 cm
- Size Small: Mouth corner distance \approx 6 cm
- Size Regular: Mouth corner distance \approx 7 cm

The orientation scale with the respective guide marks imprinted on the foil bag can be used as an aid in the selection of the correct size. It is important not to use a rigid orientation aid for the determination of the mouth corner distance, but to follow the slightly curved course of the lips. In case of doubt, we recommend using the larger size in order to ensure a better fit. The size (R for Regular, S for Small or J for Junior) is indicated on the packaging as well as on the left tab of the OptraGate. In order to ensure optimum function, we recommend inserting the largest possible size that fits the patient's mouth.

2. OptraGate may be inserted from both the 12 o'clock position with the patient lying down, or from the 8 o'clock position with the patient sitting upright (or, the 4 o'clock position, if the operator is left-handed).

- To ease placement, with the tabs of the OptraGate pointing downward, the operator should grasp the thicker intraoral ring by reaching through the thinner extraoral ring. The inner ring should be held between the thumb and middle finger and pressed together slightly (see Figure 2).


Fig 2: Correct hold on the intraoral ring with thumb and middle finger. The intraoral ring is slightly compressed in the process.

- The slightly compressed intraoral ring is then inserted into the buccal corridor on one side, so that the elastic component embraces the corner of the mouth and the outer ring is still located extraorally (see Figure 3).


Fig 3: The thicker intraoral ring is positioned between the teeth and the corner of the mouth.

- As soon as OptraGate is securely in place on one side, the other side is inserted in the same manner by slightly bending it (see Figure 4).

Fig 4: The intraoral ring is positioned behind the two corners of the mouth.

- Subsequently, the intraoral ring is placed behind the lower and upper lips and thus OptraGate assumes its final, stable position (see Figure 5). The intraoral ring can be inserted with more ease if the patient's mouth is relaxed.

Fig 5: Final positioning behind the lower and upper lip.

If, in individual cases, the intraoral ring should show a tendency to slip out of the gingivobuccal fold when the mouth is closed completely, positioning of the intraoral ring deeper inside the vestibule is usually sufficient to solve the problem. Sometimes, the use of another size may also be a solution. The securely positioned OptraGate improves the overall view and accessibility of the working field. Lateral movements of the mandible and opening of the mouth to various degrees may provide additional space and access to the treatment area. The OptraGate may even be left in place when checking the occlusion.

- To remove OptraGate, grasp the upper area of the extraoral ring with the help of a paper towel and pull it slightly downwards so that the intraoral ring is loosened from the upper gingivobuccal fold (Figure 6). OptraGate can then be easily removed from the lower gingivobuccal fold and hygienically disposed of using the paper towel.

Fig 6: Removal of the OptraGate

3 Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclar.com, and your responsible competent authority.
- The current Instructions for Use and explanation of the symbols are available on the website: www.ivoclar.com.

Warnings

- Observe the Safety Data Sheet (SDS) (available at www.ivoclar.com).
- During bleaching procedures, care should be taken to prevent any bleaching agent from getting trapped between the OptraGate and the soft tissue in order to prevent irritation of the skin and mucous membrane.

Disposal information

Remaining stocks must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks.

The following clinical residual risks are known:

- Pressure sores

4 Shelf life and storage

- Storage temperature 2–28 °C (36–82 °F)
- Keep away from sunlight.
- Do not use the product after the indicated expiry date.
- Expiry date: See note on packaging.
- Before use, visually inspect the packaging and the product for damage. In case of any doubts, please contact Ivoclar Vivadent AG or your local dealer.

5 Additional information

Keep material out of the reach of children!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the materials for their suitability and use for any purpose not explicitly stated in the Instructions.

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