Te-Econom Line

Direct restorative system





Te-Econom Line

The components of the Te-Econom line of direct restorative products are ideally coordinated and attractively priced. They include the following:

Eco-Etch® Te-Econom Flow® Te-Econom Bond® Te-Econom Plus®

Te-Econom Plus

The sculptable universal composite has been successfully used for many years. It has outstanding qualities and it is exceptionally economical to use.

Excellent physical properties

The special combination of glass fillers and spherical mixed oxides ensures high strength and low wear of the restorations.

High radiopacity

The radiopaque ytterbium trifluoride and barium glass fillers allow Te-Econom Plus to be clearly distinguished from the natural tooth structure.

Medium translucency of all shades

The translucency of the material is between that of dentin and enamel. As as a result, one material is sufficient for replacing damaged tooth structure. The shade range comprises two flowable materials in shades A2 and A3 as well six sculptable materials: A1, A2, A3, A3.5, B2 and C3.

Posterior tooth restoration

using shade A3 of the Te-Econom line of products




Initial situation, defective occlusal-distalpalatal restoration in tooth 16



Stained tooth 16 after preparation



Final situation at two-week recall





More than Some million restorations placed worldwide

Delivery forms

Te-Econom Plus Intro Pack – Article No. 610908AN

Te-Econom Plus A2, A3, Te-Econom Bond, Eco-Etch, various accessories

Te-Econom Plus Intro Pack – Article No. 610909AN

Te-Econom Plus A2, A3, A3.5, B2, Te-Econom Bond, Eco-Etch, various accessories

Te-Econom Plus Intro Pack – Article No. 648186AN

Te-Econom Plus A1, A2, A3, B2, Te-Econom Bond, Eco-Etch, various accessories

Te-Econom Plus System Pack – Article No. 610907AN

Te-Econom Plus A1, 3 x A2, 2 x A3, 1 x A3.5, B2, Te-Econom Bond, Eco-Etch, various accessories



www.ivoclarvivadent.com