

ADHESIVES

DenFil Etchant-37

Semi-gel Type Etching Agent with Phosphoric Acid



Semi-gel type enables precise application to targeted areas



No residue remains after rinsing



No clogging of the tip

INDICATIONS



Acid etching of enamel and dentin



Before Etching



After Etching

※ Microphotographs of premolar x 1,000



PACKAGE

All products			Product code	MFDS Reg. CODE
DenFil Etchant-37	Bulk Kit	5ml Etchant x 5Syringes, Disposable tips (Blue)	ET125-B000	DenFil Etchant-37
	Eco Kit	3ml Etchant x 5Syringes, Disposable tips (Blue)	ET115-E000	
		3ml Etchant x 12Syringes, Disposable tips (Blue)	ET136-K000	
		30ml Etchant x 1Syringe, Disposable tips (Blue)	ET130-B000	