

RESTORATIVES

DenFil Bulk Flow

Light-cured Flowable Bulk Fill Restorative Material



One increment
up to 4mm



Chameleon effect that ensures
a natural shade blend
in with surrounding color



Self
Leveling



Minimizing
polymerization
shrinkage

INDICATIONS



Posterior and anterior
restoration



Undercut
Blockout




Base /
Liner




PHYSICAL PROPERTIES


Flexural Strength (MPa)



Polymerization shrinkage (%)



Shrinkage Stress (MPa)



PACKAGE

PRODUCT CODE & MFDS Reg. CODE

[View details ▶](#)



DenFil Bulk Flow Kit	2g x 1 syringe	U, A1, A2, A3, B1
	2g x 2 syringes	U, A1, A2, A3, B1

SHADES GUIDE

Body Shade
A1, A2, A3,
B1, U



Flowable
Resin

POSTERIOR

RESTORATIVES

DenFil Bulk Flow

Light-cured Flowable Bulk Fill Restorative Material



PACKAGE__PRODUCT CODE & MFDS Reg. CODE

- **DenFil Bulk Flow** (2g x 1syringe or 2g x 2syringes)
- Shade: U, A1, A2, A3, B1
- Accessories: Disposable tips

ALL PRODUCTS			PRODUCT CODE	MFDS Reg. CODE
DenFil Bulk Flow Kit	2g x 1 syringe	U	BR102-RU0T	DenFil BulkFlow MTU
		A1	BR102-RA1T	DenFil BulkFlow MTA1
		A2	BR102-RA2T	DenFil BulkFlow MTA2
		A3	BR102-RA3T	DenFil BulkFlow MTA3
		B1	BR102-RB1T	DenFil BulkFlow MTB1
	2g x 2 syringes	U	BR104-RU0T	DenFil BulkFlow MTU
		A1	BR104-RA1T	DenFil BulkFlow MTA1
		A2	BR104-RA2T	DenFil BulkFlow MTA2
		A3	BR104-RA3T	DenFil BulkFlow MTA3
		B1	BR104-RB1T	DenFil BulkFlow MTB1