

RESTORATIVES

DenFil

Light-cured Restorative Hybrid
Composite Resin



INDICATIONS




Posterior and anterior
restoration




PHYSICAL PROPERTIES


Flexural Strength (MPa)



Compressive strength (MPa)



Polymerization shrinkage (%)



SHADES GUIDE

Body Shade

A1, A2, A3, A3.5, A4, B1,
B2, B3, C2, C3, D3

Opaque Shade

A20, A30

Special Shade

E, I, P


PACKAGE

PRODUCT CODE & MFDS Reg. CODE

View details ►



DenFil Kit	A Pack	A1, A2, A3, A3.5, B2
	B Pack	A2, A3, A3.5, B2, B3
	C Pack	A1, A2, A3, A3.5, C2
DenFil Refill	A1, A2, A3, A3.5, A4, B1, B2, B3, C2, C3, A20, A30, I, P, E, D3	
DenFil Capsule	Kit, A1, A2, A3, A3.5, A4, B1, B2, B3, C2, C3, A20, A30, I, P, E, D3	



DenFil

Light-cured Restorative Hybrid
Composite Resin



PACKAGE__PRODUCT CODE & MFDS Reg. CODE

- **Kit** (4g x 5syringes)
 - Shade: A pack (A1, A2, A3, A3.5, B2)
B pack (A2, A3, A3.5, B2, B3)
C pack (A1, A2, A3, A3.5, C2)
 - 5ml Bonding agent (BC Plus) x 1bottle
+ 3ml Etchant x 1syringe
 - Accessories: Brushes, Instrument, Mixing well, Mixing pad,
Disposable tips (Blue)
- **Refill** (4g x 1syringe)
 - Shade
: A1, A2, A3, A3.5, A4, B1, B2, B3, C2, C3, D3, A20, A30, I, P, E
- **Capsule** (0.25g x 20ea)
 - Shade
: A1, A2, A3, A3.5, A4, B1, B2, B3, C2, C3, D3, A20, A30, I, P, E

ALL PRODUCTS			PRODUCT CODE	MFDS Reg. CODE
DenFil Kit	A Pack	A1, A2, A3, A3.5, B2	DF104-A000	DenFil 2200
	B Pack	A2, A3, A3.5, B2, B3	DF104-B000	DenFil 2201
	C Pack	A1, A2, A3, A3.5, C2	DF104-C000	DenFil 2202
DenFil Refill	A1		DF104-RA11	DenFil 2100A1
	A2		DF104-RA22	DenFil 2100A2
	A3		DF104-RA33	DenFil 2100A3
	A3.5		DF104-RA35	DenFil 2100A3.5
	A4		DF104-RA44	DenFil 2100A4
	B1		DF104-RB11	DenFil 2100B1
	B2		DF104-RB22	DenFil 2100B2
	B3		DF104-RB33	DenFil 2100B3
	C2		DF104-RC22	DenFil 2100C2
	C3		DF104-RC33	DenFil 2100C3
	A20		DF104-ROA2	DenFil 2100A20
	A30		DF104-ROA3	DenFil 2100A30
	I		DF104-RI00	DenFil 2100I
	P		DF104-RP00	DenFil 2100P
	E		DF104-RE00	DenFil 2100E
	D3		DF104-RD33	DenFil 2100D3
DenFil Capsule	Kit		DF205-A000	-
	A1		DF205-RA11	DenFil 2101A1
	A2		DF205-RA22	DenFil 2101A2
	A3		DF205-RA33	DenFil 2101A3
	A3.5		DF205-RA35	DenFil 2101A3.5
	A4		DF205-RA44	DenFil 2101A4
	B1		DF205-RB11	DenFil 2101B1
	B2		DF205-RB22	DenFil 2101B2
	B3		DF205-RB33	DenFil 2101B3
	C2		DF205-RC22	DenFil 2101C2
	C3		DF205-RC33	DenFil 2101C3
	A20		DF205-ROA2	DenFil 2101A20
	A30		DF205-ROA3	DenFil 2101A30
	I		DF205-RI00	DenFil 2101I
	P		DF205-RP00	DenFil 2101P
	E		DF205-RE00	DenFil 2101E
	D3		DF205-RD33	DenFil 2101D3