

ADHESIVES

U-Bond Ortho

Light-cured Orthodontic Adhesive



Excellent workability with color change before and after light curing



Easy removal of remaining resin during de-bonding



Excellent adhesion to ceramic and metal brackets

INDICATIONS




Direct bonding of ceramic and metal brackets



PHYSICAL PROPERTIES

Shear bond strength (MPa)



Ensures stable bracket bonding with optimal bond strength.



Excellent workability



Before light curing



After light curing

[Patent Number : KR100932889]

The color change after light curing makes it easy to remove excess resin during bracket bonding.

PACKAGE

Type	All products	Product code	MFDS Reg. CODE
Kit	4g x 2syringes, 6ml Primer x 1bottle, Instrument, Brushes	U0108-SKK0	U0108-SKK0
		U0108-CLK0	U0108-CLK0
Refill	4g x 1syringe	U0104-SKR0	U-Bond Ortho NS
		U0104-CLR0	U-Bond Ortho NC
Primer	6ml Primer x 1bottle	OP406-R000	U-Bond Primer
	3ml Primer x 1bottle	OP403-R000	U-Bond Primer Medium