



Vonflex S Bite

Thixotropic Bite Registration Material



High tear strength



Excellent detail reproduction




Bubble-free



INDICATIONS

- Bite registration material

PHYSICAL PROPERTIES



▲ Vonflex S Bite provides highly precise and detailed occlusal reproduction.

Mixing ratio (Base:Catalyst)	1:1
Working Time (23℃)	45"
Time in mouth	90"

PACKAGE

All Products		Product CODE	MFDS Reg. CODE
Vonflex S Bite	50ml x 2cartridges, w/mixing tips	VF310-SBNB	Vonflex S Bite 0202