

Impressions

Vonflex S Light

Excellent Hydrophilic Impression Material



Excellent detail reproduction



Sufficient working time



Short setting time



High tear strength



High elastic recovery

Vonflex S Light XLV




Vonflex S Light



INDICATIONS

- Taking impression for crown and bridge, inlay, onlay

PHYSICAL PROPERTIES




Sufficient working time & Short setting time



▲ Excellent wettability for precise detail reproduction

Mixing ratio (Base:Catalyst)	1:1
Working Time (23°C)	2' 30"
Time in mouth	4' 00"



▲ Excellent flowability

PACKAGE

All Products			Product CODE	MFDS Reg. CODE
Vonflex S Light	Normal	50ml x 2cartridges, w/mixing tips	VF310-SLN0	Vonflex S Light 0102
		50ml x 12cartridges	VF360-SLN0	Vonflex S Light 0112
		50ml x 48cartridges	VF324-SLN0	Vonflex S Light 0148
	Fast	50ml x 2cartridges, w/mixing tips	VF310-SLF0	Vonflex S Light Q 0102
		50ml x 12cartridges	VF360-SLF0	Vonflex S Light Q 0112
		50ml x 48cartridges	VF324-SLF0	Vonflex S Light Q 0148
Vonflex S Light XLV	Normal	50ml x 2cartridges, w/mixing tips	VF310-SLNX	Vonflex S Light XLV 1402
		50ml x 12cartridges	VF360-SLNX	Vonflex S Light XLV 1412
		50ml x 48cartridges	VF324-SLNX	Vonflex S Light XLV 1448