

Impressions

Vonflex S Medium

Proven One Step Material Based on A-silicone



Excellent hydrophilicity



Ideal viscosity




High elastic recovery

INDICATIONS

- Taking impression for crown and bridge, inlay, onlay



PHYSICAL PROPERTIES



▲ Highly detailed and accurate impressions

Mixing ratio (Base:Catalyst)	1:1
Working Time (23°C)	2' 15"
Time in mouth	4' 00"

PACKAGE

ALL PRODUCTS			PRODUCT CODE	MFDS Reg. CODE
Vonflex S Medium	Normal	50ml x 2cartridges, w/mixing tips	VF310-SMN0	Vonflex S Medium 2502