

Impressions

Vonflex S Putty

Proven And Reliable Precision Impression Material



Convenient to mix



Bubble-free



Low oil content



Dimensional stability




INDICATIONS



- For use with individual trays
- Taking impression for crown and bridge

PHYSICAL PROPERTIES

Strain in Compression (%)



Elastic recovery (%)



▲ Easy removal due to excellent elastic recovery and high deformation recovery

Mixing ratio (Base:Catalyst)	1:1
Working Time (23°C)	2' 00"
Time in mouth	4' 00"



▲ Low oil content for improved mixability

PACKAGE

All Products			Product CODE	MFDS Reg. CODE
Vonflex S Putty	Normal	280ml Catalyst x 1ea, 280ml Base x 1ea, Measuring scoop (2ea)	VF356-SPNS	Vonflex S Putty Soft 13
	Fast		VF356-SPFS	Vonflex S Putty Soft Q 13