

ENDODONTICS

Well-Paste

Water-based Calcium Hydroxide Paste
with Barium Sulfate



Radiopaque



Optimal
flowability



Excellent accessibility
to the apex



Excellent
biocompatibility

INDICATIONS



Temporary
root canal filling



Treatment for
an infected root canal




Pulp capping
of exposed pulp




PHYSICAL PROPERTIES


Flow (mm)



pH



Film Thickness (μm)



PACKAGE

ALL PRODUCTS		PRODUCT CODE	MFDS Reg. CODE
Well-Paste	2g x 1syringe, Disposable tips, Extension holder	WS102-K000	Well-Paste Type A

Radiopacity of Well-Paste

