

ENDODONTICS

Well-Root PT

Premixed Bioceramic Putty



Low solubility



Superior compressive strength and push-out bond strength



Contains excellent biocompatible materials



Antibacterial effects by alkaline pH

INDICATIONS



Pulp capping



Repair of perforation



Root-end filling




Apexification




PHYSICAL PROPERTIES


Push-out Bond Strength (MPa)



Vickers Hardness (Hv)



Compressive Strength (MPa)



PACKAGE

ALL PRODUCTS		PRODUCT CODE	MFDS Reg. CODE
Well-Root PT	0.25g x 10capsules	WT202-K000	Well-Root PT

Properties of Well-Root PT

Properties measured according to ISO 6876:2012

Film thickness (μm)

38.8 $\pm$ 4

Radiopacity (mm Al)

5.4

Setting time (min.)

Initial 5 $\pm$ 2min Final 45 $\pm$ 5min