

ENDODONTICS

Well-Root ST

Premixed Bioceramic Sealer



Low solubility



Hydroxyapatite generation



Tight sealing



Antibacterial effects by alkaline pH

INDICATIONS




Permanent obturation of the root canal




PHYSICAL PROPERTIES

Push-out bond strength (MPa)





High pH



PACKAGE

ALL PRODUCTS		PRODUCT CODE	MFDS Reg. CODE
Well-Root ST	2g x 1syr. with Disposable tip(20ea) & Extension holder	WR102-K00C	Well-Root ST Type C

Tight sealing



Optical microscope image of microleakage test (dye penetration test); No dye penetration